



Product: HUMAN FACTOR V  
Catalog #: IHFV  
Lot Number: Sample  
Concentration: 7.3 mg/mL  
Specific Activity: 44 U/mg\*  
Extinction Coefficient:  $E_{1\text{cm}, 280\text{nm}}^{1\%} = 9.6$   
Molecular Weight: 330000  
Buffer: 50% Glycerol/ H<sub>2</sub>O (vol/vol)  
Recommended Storage: -20 ° C

Amount Shipped: \_\_\_\_\_  
Date Shipped: \_\_\_\_\_  
Shipped by: \_\_\_\_\_  
Expiration: \_\_\_\_\_

**PRODUCT/CAUTION NOTES:**

Products contained within this shipment are of human blood/plasma origin. Although the starting material was tested prior to initiation of the manufacturing process, and was found negative or nonreactive for anti-HIV-1/2, HIV-1 antigen(s), HBsAg, STS, anti-HCV, anti-HBcore and anti-HTLV I & II, extreme caution should be used when handling this material as there is a margin of error in all tests.

This product is to be used for IN VITRO research purposes only, and is not to be used for clinical, diagnostic or human use!

\*Determined by Factor V clotting assay. One (1) unit is equivalent to the Factor V activity in one (1) milliliter of normal human plasma.