



Certificate of Analysis

Lot: SAMPLE SAMPLE

Description: Porcine Thrombin

Concentration: 1.11 mg/ml

Aliquots: 1 x 0.10 mg

Volume/aliquot: 0.090 ml

Activity: 2054.00 NIH u/mg

Total protein: 0.10 mg

Buffer: 50 mM sodium citrate/ 0.2 M NaCl/
0.1% PEG-8000/ pH 6.5

Storage: < -60 °C

Expiration Date: 09/2006

The Porcine Thrombin was activated from homogeneous porcine prothrombin by activation with Bovine Factor Xa, Factor Va, and phospholipid. These activating enzymes were removed after activation.