

NBP1-00085G**Beclin 1 Antibody (4H10) [DyLight 488]****Unit Size**

0.1 ml

Concentration

1.2 mg/ml

Storage

Store at 4C in the dark.

Applications

IHC-P, WB

Species

Bovine, Chicken, Equine, Human, Mouse, Porcine, Primate, Rat

Host

Mouse

Clonality

Monoclonal

Isotype

IgG1 Kappa

Conjugate

DyLight 488

Purity

IgG purified

Preservative

No Preservative

Specificity/Sensitivity

Beclin 1 (4H10)

Immunogen

Full length human Beclin 1 protein [Swiss-Prot# Q14457]

Recommended Dilutions

Western Blot 2.0 ug/ml, Immunohistochemistry-Paraffin 1:50- 1:100

Buffer

50mM sodium borate

Uses

This antibody is useful for Western blot, where a band is seen ~52 kDa and Immunohistochemistry-Paraffin.

References (more available at www.novusbio.com/NBP1-00085G)

Klionsky, D.J., Emr, S. Autophagy as a regulated pathway of cellular degradation. Science. 290: 1717-1721, 2000.

Images**This product does not currently have an image. Earn rewards by submitting your image of this product. Submit a review with your image here:****www.novusbio.com/NBP1-00085G**

Updated 2/13/2012 2.0

*This product is for research use only and is not approved for use in humans or in clinical diagnosis. Products are guaranteed for six months from the date of receipt except for peptides and proteins which are guaranteed for 3 months. For more information on our guarantee, please visit www.novusbio.com/guarantee.***Novus USA**8100 Southpark Way, A-8
Littleton, CO 80120
p) 888-506-6887
p) 303-730-1950
f) 303-730-1966
novus@novusbio.com**Novus Canada**461 North Service Road West,
Unit B37
Oakville, ON L6M 2V5
p) 855-668-8722
f) 905-827-6402
canada@novusbio.com**Novus Europe**12 Cambridge Science Park
Cambridge, UK CB4 0FQ
UK: p) +44-1223-426-001, f) +44-871-971-1635
DE: p) +49-6922-22340-60, f) +49-0800-58926-79
IT: p) +39-02-4032-6786, f) +39-02-4032-6340
FR: p) +33-1-76-77-45-30, f) +33-1-76-77-45-31
europe@novusbio.com