

DESCRIPTION

Source	<i>E. coli</i> -derived Ser19-Ser178 Accession # Q9Z1Y5
N-terminal Sequence Analysis	Ser19
Structure / Form	Noncovalently-linked homodimer
Predicted Molecular Mass	18.5 kDa (monomer)

SPECIFICATIONS

SDS-PAGE	18.5 kDa, reducing conditions
Activity	Measured in a cell proliferation assay using MC/9-2 mouse mast cells. Thompson-Snipes, L. <i>et al.</i> (1991) <i>J. Exp. Med.</i> 173 :507. The ED ₅₀ for this effect is typically 1.5-6 ng/mL.
Endotoxin Level	<1.0 EU per 1 µg of the protein by the LAL method.
Purity	>95%, by SDS-PAGE under reducing conditions and visualized by silver stain.
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with BSA as a carrier protein. See Certificate of Analysis for details.

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 100 µg/mL in sterile PBS containing at least 0.1% human or bovine serum albumin.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 3 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Interleukin 10, also known as cytokine synthesis inhibitory factor (CSIF), is the charter member of the IL-10 family of α -helical cytokines that also includes IL-19, IL-20, IL-22, and IL-24 (1, 2). IL-10 is secreted by many activated hematopoietic cell types as well as hepatic stellate cells, keratinocytes, and placental cytotrophoblasts (2 - 5). Mature guinea pig IL-10 shares 75% - 84% amino acid sequence identity with bovine, canine, equine, feline, human, mouse, ovine, porcine, and rat IL-10. IL-10 is a 178 amino acid molecule that contains two intrachain disulfide bonds and is expressed as a noncovalently associated homodimer (6 - 8). The IL-10 dimer binds to two IL-10 R α /IL-10 R1 chains, resulting in recruitment of two IL-10 R β /IL-10 R2 chains and activation of a signaling cascade involving JAK1, TYK2, and STAT3 (9). IL-10 R β does not bind IL-10 by itself but is required for signal transduction (1). IL-10 R β also associates with IL-20 R α , IL-22 R α , or IL-28 R α to form the receptor complexes for IL-22, IL-26, IL-28, and IL-29 (10 - 12). IL-10 is a critical molecule in the control of viral infections and allergic and autoimmune inflammation (13 - 15). It promotes phagocytic uptake and Th2 responses, but suppresses antigen presentation and Th1 proinflammatory responses (2).

References:

1. Pestka, S. *et al.* (2004) *Annu. Rev. Immunol.* **22**:929.
2. O'Garra, A. and P. Vieira (2007) *Nat. Rev. Immunol.* **7**:425.
3. Mathurin, P. *et al.* (2002) *Am. J. Physiol. Gastrointest. Liver Physiol.* **282**:G981.
4. Grewe, M. *et al.* (1995) *J. Invest. Dermatol.* **104**:3.
5. Szony, B.J. *et al.* (1999) *Mol. Hum. Reprod.* **5**:1059.
6. Scarozza, A.M. *et al.* (1998) *Cytokine* **10**:851.
7. Windsor, W.T. *et al.* (1993) *Biochemistry* **32**:8807.
8. Syto, R. *et al.* (1998) *Biochemistry* **37**:16943.
9. Kotenko, S.V. *et al.* (1997) *EMBO J.* **16**:5894.
10. Kotenko, S.V. *et al.* (2000) *J. Biol. Chem.* **276**:2725.
11. Hor, S. *et al.* (2004) *J. Biol. Chem.* **279**:33343.
12. Sheppard, P. *et al.* (2003) *Nat. Immunol.* **4**:63.
13. Fitzgerald, D.C. *et al.* (2007) *Nat. Immunol.* **8**:1372.
14. Wu, K. *et al.* (2007) *Cell. Mol. Immunol.* **4**:269.
15. Blackburn, S.D. and E.J. Wherry (2007) *Trends Microbiol.* **15**:143.