

DESCRIPTION

Species Reactivity	Human
Specificity	Detects human Complement Factor H-related 1/CFHR1 in direct ELISAs and Western blots. In direct ELISAs and Western blots, no cross-reactivity with recombinant human (rh) CFHR2 or rhCFHR5 is observed.
Source	Monoclonal Mouse IgG _{2B} Clone # 442127
Purification	Protein A or G purified from hybridoma culture supernatant
Immunogen	Mouse myeloma cell line NS0-derived recombinant human Complement Factor H-related 1/CFHR1 Glu19-Ala328 Accession # Q03591
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details.

APPLICATIONS

Please Note: Optimal dilutions should be determined by each laboratory for each application. *General Protocols* are available in the *Technical Information* section on our website.

	Recommended Concentration	Sample
Western Blot	1 µg/mL	Recombinant Human Complement Factor H-related 1/CFHR1 under non-reducing conditions
Immunoprecipitation	25 µg/mL	Conditioned cell culture medium spiked with Recombinant Human Complement Factor H-related 1/CFHR1, see our available Western blot detection antibodies

PREPARATION AND STORAGE

Reconstitution	Reconstitute at 0.5 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below.
Stability & Storage	<p>Use a manual defrost freezer and avoid repeated freeze-thaw cycles.</p> <ul style="list-style-type: none"> ● 12 months from date of receipt, -20 to -70 °C as supplied. ● 1 month, 2 to 8 °C under sterile conditions after reconstitution. ● 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

Complement factor H-related protein 1 (CFHR1) is a 43 kDa, secreted member of the factor H family of glycoproteins. It is produced by hepatocytes and circulates as two differentially glycosylated isoforms (37 kDa and 43 kDa). Mature human CFHR1 is 312 aa in length. It contains five approximately 60 aa SCRs (short consensus repeats/CCPs/SUSHI repeats) that basically constitute the entire molecule. CFHR1 may play a role in lipoprotein complexes that bind LPS to neutrophils. There is no reported rodent counterpart to CFHR1. Over aa 19-143 of the CFHR1 precursor, human CFHR2 and CFHR5 show 98% and 86% aa identity, respectively.