


Polyclonal Antibody to Serotonin (5-HT)

Alternate names: 5 HT, 5 hydroxytryptamine
Catalog No.: 20080
Quantity: 0.1 ml
Host: Rabbit
Pictures:



Staining with 5-HT (Serotonin) Rabbit Antibody 20080

For research and in vitro use only. Not for diagnostic or therapeutic work.

Material Safety Datasheets are available at www.acris-antibodies.com or on request.

Antibody Hotline - Technical Questions - Antibody Location Service
Free Call: 0800-2274746 (Germany only) - www.acris-antibodies.com