

This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännas | Sweden | +46 035 33033 | www.agrisera.com

HTR12 | histone H3-like centromeric protein (200 µl)

400 •

AS12 2226-200 | clonality: **polyclonal** | host: **rabbit** | reactivity: ***Arabidopsis thaliana***

[PRODUCT INFORMATION IN PDF](#)

Qty:

[Add to cart](#)

[Add to wishlist](#)

[Add review](#)

Availability:

0 st

Delivery:

3-6 working days

Shipping:

World-wide

Item No:

AS12 2226-200

Info: [More information](#) [Product suggestions](#) [Add review](#)

[Plant protein extraction buffer](#)

[Secondary antibodies](#)