

## Anti-Fd-GOGAT (Ferredoxin-dependent glutamate synthase, chloroplastic) antibody, rabbit polyclonal

81-025 100 μg

**Storage:** Ship at 4°C and store at -20°C. Do not freeze.

**Immunogen:** Purified recombinant maize Fd-GOGAT protein, full-size, no-tag attached **Reactivity:** Fd-GOGAT including those of maize, arabidopsis, spinach, and cyanobacterium (Synechocystis)

**Validation:** Specificity has been validated by WB with recombinant full-size maize Fd-GOGAT protein

## Applications:

- 1. Western blotting (1/2,000-1/5,000 dilution)
- 2. ELISA (assay dependent)

Other applications have not been tested.

**Purity:** IgG, purified with protein A/G mix.

**Form:** 2 mg/ml in PBS, 50% glycerol. Filter sterilized. No preservative or carrier protein added.

Background: Glutamine oxoglutarate aminotransferase (abbreviated as GOGAT) is an enzyme involved in synthesis of glutamate from glutamine and alpha-ketoglutarate. GOGAT has two forms in plants: ferredoxin-dependent GOGAT (Fd-GOGAT) and NADH-dependent GOGAT (NADH-GOGAT). 95% of GOGAT found in plants is the Fd-GOGAT type. Fd-GOGAT is encoded by two genes, glu1 and glu2 in Arabidopsis. Fd-GOGAT (both forms) is highly conserved among plants, red algae, and cyanobacteria. Ferredoxin-dependent glutamate synthase, chloroplastic (Fd-GOGAT) is involved in glutamate biosynthesis in leaf. This protein required for the reassimilation of ammonium ions generated during photorespiration. Gene name is GlsF.

**Data Link**: UniProtKB: <u>P23225</u> (Z. mays), <u>Q51579</u> (P. boryanum), <u>P55038</u> (*Synechocystis sp.*)






Fig.1 Western Blot of Fd-GOGAT in plant leaf extract.

- 1. Full-size recombinant maize Fd-GOGAT protein
- 2. Arabidopsis leaf extract,  $10 \mu g$
- 3. Maize leaf extract, 10 µg

Anti-Fd-GOGAT antibody was used at 1/2,500 dilution. Second antibody (goat anti-rabbit IgG antibody HRP-conjugated, ab97051) was used at 1/10,000 dilution.

Molecular masses of maize and arabidopsis Fd-GOGAT are 175 kDa and 168 kDa, respectively..




Fig.2 Western Blot of Fd-GOGAT in Cyanobacterium, Synechocystis PCC 6803

Sample: Crude extract of The Fd-GOGAT antibody was used at 1/2000 dilutions Molecular mass is 169 kDa

**Reference**: This antibody has been used in the following publications.

- 1 Sakakibara H. et al. Molecular cloning and characterization of complementary DNA encoding for ferredoxin-dependent glutamate synthase in maize leaf. J Biol Chem. 1991 Feb 5;266(4):2028-35. PMID: 1989968. WB; maize
- 2. Sakakibara H., Kawabata S., Hase T. and Sugiyama T. (1992) Differential effects of nitrate and light on the expression of glutamine synthetase and ferredoxin-dependent glutamate synthase in maize. Plant Cell Physiol., 33, 1193-1198. Googl Scholar: abstract/33/8/1193/1860644 WB; maize
- 3.Kimata-Ariga Y and Hase T. Multiple complexes of nitrogen assimilatory enzymes in spinach chloroplasts: possible mechanisms for the regulation of enzyme function. PLoS One. 2014 Oct 1;9(10):e108965. PMID: 25271437 WB;spinach