


TubeSpin® Bioreactor 50 - 180

DESCRIPTION

50 ml tubes with membrane filter TubeSpin® Bioreactor 50, TPP® 180 pcs.

SPECIFICATIONS

Number of tubes in set	180
Volume	50 ml



CAT. NUMBER

BS-010158-CK 50 ml tubes with membrane filter
 TubeSpin® Bioreactor 50, TPP®
 180 pcs.

ACCESSORIES



RTS-1
Including TPP TubeSpin® Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1 is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)



RTS-1C
Including TPP TubeSpin® Bioreactor vessels 50ml, 20pcs
Personal bioreactor

RTS-1C is personal bioreactor which utilize patented Reverse-Spin® technology that applies non-invasive, mechanically driven, low energy consumption, innovative type of agitation where cell suspension is mixed by the singleuse falcon ...

[read more](#)