



The Immunogold Specialist

Product Analysis

EM Goat anti-Mouse IgG+ IgM (H&L): 10nm Gold

Product Information

Product Code:	EM.GAF10
Batch Number:	****
Protein Absorption:	Minimum cross reaction to Human, Bovine and Horse serum proteins
Protein Concentration:	23.4 µg / ml
Storage Buffer:	Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN ₃ (preservative) and 20% Glycerol
Storage:	Store at 4°C

Quality Control

Protein Detected:	7 pg after silver enhancement	Passed
Optical Density:	3.1	Passed
Mean Diameter:	9.6nm	Passed
Coefficient of Variation:	7.8%	Passed
Percentage Single Particles:	96%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBinternational

Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: info@bbigold.com

www.bb-international.com

A complete catalogue of products can be found on our website or CD Rom