



The Immunogold Specialist

## **Product Analysis**

### **EM Goat anti-Rabbit IgG (H&L): 10nm Gold**

#### Product Information

Product Code: EM.GAR10  
Batch Number: \*\*\*\*  
Protein Absorption: Minimum cross reaction to Human serum proteins  
Protein Concentration: 23.2 µg / ml  
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA, 15mM NaN<sub>3</sub> (preservative) and 20% Glycerol  
Storage: Store at 4°C

#### Quality Control

Protein Detected:	2pg after silver enhancement	<b>Passed</b>
Optical Density:	3.1	<b>Passed</b>
Mean Diameter:	9.4nm	<b>Passed</b>
Coefficient of Variation:	6.0%	<b>Passed</b>
Percentage Single Particles:	98%	<b>Passed</b>

QC approved: \_\_\_\_\_

**For in vitro research only. Not for diagnostic use**

**To reorder or for further technical support:**

**BBInternational**  
Golden Gate, Ty Glas Avenue,  
Llanishen, Cardiff, CF14 5DX, UK  
Tel: 44 (0) 2920 747 232  
Fax: 44 (0) 2920 747 242  
E-mail: [info@bbigold.com](mailto:info@bbigold.com)  
[www.bb-international.com](http://www.bb-international.com)

***A complete catalogue of products can be found on our website or CD Rom***