



The Immunogold Specialist

## **Product Analysis**

### **EM Goat anti-Rat IgG (H&L): 10nm Gold**

#### Product Information

Product Code: EM.GAT 10

Batch Number: \*\*\*\*

Protein Absorption: Minimum cross reaction to Human, Bovine and Horse serum proteins

Protein Concentration: 23.2 µg/ml

Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA. 15mM NaN<sub>3</sub> (preservative) and 20% Glycerol

Storage: Store at 4°C

#### Quality Control

Protein Detected:	2pg after silver enhancement	<b>Passed</b>
Optical Density:	3.0	<b>Passed</b>
Mean Diameter:	10.8nm	<b>Passed</b>
Coefficient of Variation:	11.2%	<b>Passed</b>
Percentage Single Particles:	98%	<b>Passed</b>

QC approved: \_\_\_\_\_

**For in vitro research only. Not for diagnostic use**

**To reorder or for further technical support:**

**BBinternational**

Golden Gate, Ty Glas Avenue,  
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: [info@bbigold.com](mailto:info@bbigold.com)

[www.bb-international.com](http://www.bb-international.com)

***A complete catalogue of products can be found on our website or CD Rom***