



The Immunogold Specialist

Product Analysis

EM Protein A: 10nm Gold

Product Information

Product Code: EM.PAG 10
Batch Number: ****
Protein Absorption: Not Applicable
Protein Concentration: 4.8µg / ml
Storage Buffer: Tris Buffered Saline pH8.2 containing 1% BSA,
15mM NaN₃ (preservative) and 20% Glycerol
Storage: Store at 4°C

Quality Control

Protein Detected:	4pg	Passed
Optical Density:	3.1	Passed
Mean Diameter:	10.1nm	Passed
Coefficient of Variation:	6.3%	Passed
Percentage Single Particles:	100%	Passed

QC approved: _____

For in vitro research only. Not for diagnostic use

To reorder or for further technical support:

BBinternational

Golden Gate, Ty Glas Avenue,
Llanishen, Cardiff, CF14 5DX, UK

Tel: 44 (0) 2920 747 232

Fax: 44 (0) 2920 747 242

E-mail: info@bbigold.com

www.bb-international.com

A complete catalogue of products can be found on our website or CD Rom