

LoTox nMus m 1

## LoTox™ Natural Mus m 1

Product Code: LTN-MM1-1

Allergen: Natural Mus m 1

Lot No: XXXXX

Source: Mouse Urine

Mol. Wt: 18-20 kD

Purification: Purified from concentrated, male, mouse urine by gel

filtration and ion exchange HPLC. Purity > 95% by

silver-stained SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative-free and carrier-free in endotoxin free

20 mM Tris-HCl, pH 8.1. Filtered through 0.22µ filter.

Storage: Store at -20°C

Notes: (1) A LoTox™ product, Endotoxin < 0.01 EU/µg.

(2) Avoid repeated Freeze/Thaw cycles.

## Allergens are provided for research and commercial use in vitro. Not for human in vivo or therapeutic use.

## **REFERENCES:**

- 1. Ferrari E, Lodi T, Sorbi RT, Tirindelli R, Cavaggioni A, Spisni A. Expression of a lipocalin in Pichia pastoris: secretion, purification and binding activity of a recombinant mouse major urinary protein. FEBS Lett. 1997 Jan 13;401(1):73-7.
- 2. Lorusso, JR, Moffat S, Ohman JL Jr. Immunological and biochemical properties of the major mouse urinary allergen (Mus m 1). J Allergy Clin Immunol 1986 Nov;78(5 Pt 1):928-37.