



Recombinant Can f 6

Product Code: RE-CF6-1

Allergen: rCan f 6 (*Canis familiaris* allergen 6.0101)

Lot No: xxxx

Source: *E. coli*

Mol. Wt: 21 kD

Purification: Multi-step chromatographic purification.
Purity > 95% by SDS-PAGE.

Concentration: See product insert

Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4.
Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert

Notes: 1. rCan f 6 has a C-Terminus 6XHis-tag.

Gel Image:
rCan f 6



For research and commercial use in vitro: not for human in vivo or therapeutic use.

An InBio® product. Made in the USA

Indoor Biotechnologies, Inc.
700 Harris Street
Charlottesville, VA 22903, USA
www.inbio.com

InBio® US: mail@inbio.com
InBio® UK: info@indoorbiotech.co.uk