



Recombinant Tri a 14

Product Code: RE-TA14-1

Allergen: rTri a 14 (*Triticum aestivum* allergen 14)

Lot No: XXXXX

Source: *E. coli*

Mol. Wt: 10.5 kD

Purification: Purified from *E. coli* culture by multi-step HPLC chromatography. Purity > 95 % by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative and carrier-free in 20mM Tris, 0.3M NaCl, 0.2M L-arginine pH 8.0. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) rTri a 14 has a C-terminal 6xHis-tag.

Gel Image:
rTri a 14



For research and commercial use in vitro: not for human in vivo or therapeutic use.

An InBio® product. Made in the USA