

rAsp f 1

# Recombinant Asp f 1.0101

# Product Code: RP-AF1-1

Allergen:	rAsp f 1 (Aspergillus fumigatus allergen 1.0101)
Lot No:	XXXXX
Source:	Pichia pastoris
Mol. Wt:	18 kD
Purification:	Purified by multi-step HPLC. Purity on silver- stained SDS-PAGE >95%.
Concentration:	See product insert.
Formulation:	Preservative-free and carrier-free in phosphate buffered saline, pH 7.4. Filtered through $0.22\mu$ filter.
Storage:	Store at -20°C. Avoid repeated Freeze/Thaw cycles.
Notes:	<ol> <li>rAsp f 1 has a C-terminal 6x His-tag.</li> <li>Upper band is Asp f 1 dimer.</li> </ol>

# An InBio<sup>™</sup> product

# For Research Use Only: Not for Diagnostic or Therapeutic Use

# **REFERENCES:**

- 1. Arruda, LK, Platts-Mills, TAE, Fox, JW and Chapman, MD (1990). Aspergillus fumigatus allergen I, a major IgE binding protein, is a member of the mitogillin family of cytotoxins. J. Exp. Med., 172:1529-1532.
- Arruda LK, Platts-Mills TAE, Longbottom JL, Chapman MD (1992). Aspergillus fumigatus: Identification of 16kd, 18kd and 45kd antigens recognized by human IgG and IgE antibodies and murine monoclonal antibodies. J. Allergy Clin. Immunol. 89:1166-1176.
- 3. Arruda LK, Mann BJ, Chapman MD (1992). Selective expression of a major allergen and cytotoxin, Asp f I, in Aspergillus fumigatus. J. Immunol.149:3354-3359.
- 4. Sporik RB, Arruda LK, Woodfolk J, Chapman MD, Platts-Mills TAE. (1993). Environmental exposure to Aspergillus fumigatus allergen 1 (Asp f I). Clin Exp Allergy 23:326-331.
- 5. EI-Dahr JM, Fink R, Selden R, Arruda LK, Platts-Mills TAE, Heymann PW. (1994). Development of immune responses to Aspergillus at an early age in children with cystic fibrosis. Am J Respir Crit Care Med 150:1513-1518.