



Recombinant Jug r 3

Product Code: RPI-JR3-1

Allergen: rJug r 3 (*Juglans regia* allergen 3.0101)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 10 kD

Gel Image:

rJug r 3



Purification: Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

Concentration: See product insert.

Formulation: Preservative and carrier-free in phosphate buffered saline, 300mM NaCl, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) rJug r 3 has a N-terminal 6xHis-tag.

For research and commercial use in vitro: not for human in vivo or therapeutic use.

An InBio® product. Made in the USA

Indoor Biotechnologies, Inc.
700 Harris Street
Charlottesville, VA 22903, USA
www.inbio.com

InBio® US: mail@inbio.com
InBio® UK: info@indoorbiotech.co.uk
InBio® India: info@inbioindia.com