



## Recombinant Mal d 1

**Product Code:** RPI-MD1-1

**Gel Image:**

**Allergen:** rMal d 1 (*Malus domestica* allergen 1.0108)

**Lot No:** XXXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 17.5 kD

**Purification:** Purified from *Pichia pastoris* culture by multi-step chromatography. Purity > 95 % by SDS-PAGE

**Concentration:** See product insert.

**Formulation:** Preservative and carrier-free in 20mM tris, 0.3 M NaCl, pH 8.0. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert.

**Notes:** (1) rMal d 1 has a N-terminal 6xHis-tag



**For research and commercial use in vitro: not for human in vivo or therapeutic use.**

**An InBio® product. Made in the USA**

Indoor Biotechnologies, Inc.  
700 Harris Street  
Charlottesville, VA 22903, USA  
www.inbio.com

InBio® US: mail@inbio.com  
InBio® UK: info@indoorbiotech.co.uk  
InBio® India: info@inbioindia.com