



Recombinant Pen a 1

Product Code: RPI-PA1-1

Allergen: rPen a 1 (*Penaeus aztecus* allergen 1.0101)

Lot No: XXXXX

Source: *Pichia pastoris*

Mol. Wt: 33 kD

Purification: Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

Concentration: See product insert.


Formulation: Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

Storage: Store at -20°C. Avoid repeated freeze-thaw cycles.

Expiration: See product insert.

Notes: (1) rPen a 1 has a C-terminal 6xHis-tag.

Gel Image:
rPen a 1



For research and commercial use in vitro: not for human in vivo or therapeutic use.

An InBio® product. Made in the USA

Indoor Biotechnologies, Inc.
700 Harris Street
Charlottesville, VA 22903, USA
www.inbio.com

InBio® US: mail@inbio.com
InBio® UK: info@indoorbiotech.co.uk
InBio® India: info@inbioindia.com