



## Recombinant Pla I 1

**Product Code:** RPI-PL1-1

**Allergen:** rPla I 1 (*Plantago lanceolata* allergen 1.0101)

**Lot No:** XXXXX

**Source:** *Pichia pastoris*

**Mol. Wt:** 15 kD

**Purification:** Purified from *Pichia* culture supernatant using multi-step chromatography. Purity >95% by SDS-PAGE.

**Concentration:** See product insert.

**Formulation:** Preservative and carrier-free in phosphate buffered saline, pH 7.4. Filtered through 0.22µm filter.

**Storage:** Store at -20°C. Avoid repeated freeze-thaw cycles.

**Expiration:** See product insert.

**Notes:** (1) rPla I 1 has a N-terminal 6xHis-tag.

**Gel Image:**  
rPla I 1



**For research and commercial use in vitro: not for human in vivo or therapeutic use.**

**An InBio® product. Made in the USA**

Indoor Biotechnologies, Inc.  
700 Harris Street  
Charlottesville, VA 22903, USA  
www.inbio.com

InBio® US: mail@inbio.com  
InBio® UK: info@indoorbiotech.co.uk  
InBio® India: info@inbioindia.com