

Product Data Sheet

Product:	PepMix Collection EBV (25ug/peptide)
Product code:	PM-C-EBV-2
Source:	Epstein-Barr Virus (EBV) Selected proteins
Protein ID:	see datasheet
No. of peptides:	2096 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	diverse
Sequence:	14 different PepMixes, 1 aliquots of 14 peptide pools each containing 15µg per peptide /tube
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months