

Product Data Sheet

Product:	PepMix Collection Prostate Cancer
Product code:	PM-C-PC
Source:	Human Selected proteins
Protein ID:	see datasheet
No. of peptides:	793 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	diverse
Sequence:	7 different TAA PepMixes + 2 Controls (15mers, overlap 11); Human (CEA), Human (Metalloreductase STEAP1), Human (Myc protein), Human (PAP), Human (Prostate-specific antigen PSA), Human (TCRgamma alternate reading frame protein TARP), Human (WT1/WT33), CEF Pool (extended), Human (Actin)
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C.
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months