

Product Data Sheet

Product:	CEFX Ultra SuperStim Pool MHC-I Subset
Product code:	PM-CEFX-4
Source:	Other/None Selected proteins
No. of peptides:	80 peptides - This pool consists of 80 known peptide epitopes for a broad range of HLA class 1 alleles and different infectious agents for T-cell stimulation of populations with a diverse ethnic background (see annotations).
Amount:	15 nmol (approx. 25µg) / peptide for stimulation of $2,5 \times 10^8$ cells
Length:	diverse
Sequence:	see additional product documentation
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months