

Product Data Sheet

Product:	PepMix Human (CTLA4)
Product code:	PM-CTLA4
Source:	Human Cytotoxic T-lymphocyte protein
Protein ID:	P16410
Gene ID:	CTLA4
No. of peptides:	53 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	223 aa
Sequence:	MACLGFQRHK AQLNLARTW PCTLLFFLLF IPVFCKAMHV AQPAVVLASS RGIASFVCEY ASPGKATEVR VTVLRQADSQ VTEVCAATYM MGNELTFLDD SICTGTSSGN QVNLTIQGLR AMDTGLYICK VELMYPPPPY LGIGNGTQIY VIDPEPCPDS DFLLWILAAV SSGLEFFYSFL LTAVSLSKML KKRSPLETTGV YVKMPPEPE CEKQFPYFI PIN
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months