

Product Data Sheet

Product:	PepMix EBV (BARF1)
Product code:	PM-EBV-BARF1
Source:	Epstein-Barr Virus (EBV) strain B95-8 Secreted protein BARF1 (p53)
Protein ID:	P03228
Gene ID:	BARF1
No. of peptides:	53 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	221 aa
Sequence:	MARFIAQLLL LASCVAAGQA VTAFLGERVT LTSYWRRVSL GPEIEVSWFK LGPGEEQVLI GRMHHDVIFI EWPFRRGFFDI HRSANTFFLV VTAANISHDG NYLCRMKLGE TEVTKQEHLIS VVKPLTILSVH SERSQFPDFS VLTVTCTVNA FPHPHVQWLM PEGVEPAPTA ANGGVMKEKD GSLSVAVDLS LPKPWHLPTV CVGKNDKEEA HGVYVSGYLS Q
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months