

Product Data Sheet

Product:	PepMix HIV-1 (NEF) Ultra
Product code:	PM-HIV-NEF
Source:	Human Immunodeficiency Virus (HIV) Nef protein (NEF)
No. of peptides:	150 peptides (Selected immunogenic peptides (15mers))
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	diverse
Sequence:	see additional product documentation
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months