

Product Data Sheet

Product:	PepMix HRSVB (Major surface glycoprotein G)
Product code:	PM-HRSVB-MSG
Source:	Respiratory Syncytial Virus (RSV) type B (strain B1) Major surface glycoprotein G
Protein ID:	O36633
Gene ID:	G
No. of peptides:	72 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	299 aa
Sequence:	MSKHKNQRTA RTLEKTWDTL NHLIVISSCL YRLNLKSIAQ IALSVLAMII STSLIIAAII FIISANHKVT LTTVTVQTIK NHTEKNITTY LTQVPPERVS SSKQPTTSP IHTNSATTSP NTKSETHHTT AQTGRTTTS TQTNKPSTKP RLKNPPKKPK DDYHFEVFNF VPCSIGNNQ LCKSICKTIP SNKPKKKPTI KPTNKPTTKT TNKRDPKTPA KTTKKETTTN PTKKPTLTTT ERDTSTSQST VLDTTTTLEHT IQQQSLHSTT PENTPNSTQT PTASEPSTSN STQNTQSHA
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months