

Product Data Sheet

Product:	PepMix HRSVB (Nucleocapsid Protein N)
Product code:	PM-HRSVB-NCPN
Source:	Respiratory Syncytial Virus (RSV) type B (strain B1) Nucleoprotein
Protein ID:	O42053
Gene ID:	N
No. of peptides:	95 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	391 aa
Sequence:	MALSKVKLND TLNKDQLLSS SKYTIQRSTG DNIDTPNYDV QKHLNKLCGM LLITEDANHK FTGLIGMLYA MSRLGREDTI KILKDAGYHV KANGVDITTY RQDINGKEMK FEVLTLSSLT SEIQVNIEIE SRKSYKKMLK EMGEVAPEYR HDSPDCGMII LCIAALVITK LAAGDRSGLT AVIRRANNVL KNEIKRYKGL IPKDIANSFY EVFEKHPHLI DVFVHFGIAQ SSTRGGSRVE GIFAGLFMNA YGSGQVMLRW GVLAKSVKNI MLGHASVQAE MEQVVEVY EY AQKLGGEAGF YHILNPKAS LLSLTQFPNF SSVVLGNAAG LGIMGEYRGT PRNQDLYDAA KAYAEQLKEN GVINYSVLDL TAEELEAIKN QLNPKEDDVE L
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months