

Product Data Sheet

Product:	PepMix Mouse (Survivin)
Product code:	PM-MM-Survivin
Source:	Mouse Baculoviral IAP repeat-containing protein 5 (Survivin)
Protein ID:	O70201
Gene ID:	Birc5
No. of peptides:	33 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to $2,5 \times 10^8$ cells)
Length:	140 aa
Sequence:	MGAPALPQIW QLYLKNYRIA TFKNWPFLD CACTPERMAE AGFIHCPTEN EPDLAQCFFC FKELEGWEPD DNPIEEHRKH SPGCAFLTVK KQMEELTVSE FLKLDQRQAK NKIAKETNNK QKEFEETAKT TRQSIEQLAA
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months