

## Product Data Sheet

<b>Product:</b>	PepMix Human (Melan-A/MART-1)
<b>Product code:</b>	PM-MelA
<b>Source:</b>	Human Melanoma antigen recognized by T-cells 1 (MART-1)
<b>Protein ID:</b>	Q16655
<b>Gene ID:</b>	MLANA
<b>No. of peptides:</b>	27 peptides
<b>Amount:</b>	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 <sup>8</sup> cells)
<b>Length:</b>	118 aa
<b>Sequence:</b>	MPREDAHFIY GYPKKGHGHS YTTAEEAAGI GILTVILGVL LLIGCWYCRR RNGYRALMDK SLHVGTQCAL TRRCPQEGFD HRDSKVSLQE KNCEPVPVNA PPAYEKLSAE QSPPPYSP
<b>Purity:</b>	Development Grade: each peptide purified to > 70% (HPLC/MS)
<b>Appearance:</b>	white to off-white lyophilisate
<b>Remark:</b>	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
<b>Storage conditions:</b>	Store at - 20°C
<b>Application protocol:</b>	<a href="https://jpt.com/support-contact/resources/application-protocols/">https://jpt.com/support-contact/resources/application-protocols/</a>

**Use recommended within 6 months**