

Product Data Sheet

Product:	PepMix Human (PD-L1)
Product code:	PM-PD-L1
Source:	Human Programmed cell death 1 ligand
Protein ID:	Q9NZQ7
Gene ID:	CD274
No. of peptides:	70 peptides
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	290 aa
Sequence:	MRIFAVFIFM TYWHLNNAFT VTVPKDLYVV EYGSNMTIEC KFPVEKQLDL AALIVYWEME DKNIIQFVHG EEDLKVQHSS YRQRARLLKD QLSLGNAALQ ITDVKLQDAG VYRCMISYGG ADYKRITVKV NAPYNKINQR ILVVDPTSE HELTCQAEGY PKAEVIWTSS DHQVLSGKTT TTNSKREEKL FNVSTLRIN TTTNEIFYCT FRRLDPEENH TAEIPELPEL LAHPPNERTH LVILGAILLC LGVALTFIFR LRKGRMMDVK KCGIQDTNSK KQSDTHLEET
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months