

Product Data Sheet

Product: PepMix P. vivax (CSP) Ultra

Label name: P. vivax (CSP) Ultra

Product code: PM-PLVI-CSP-1

Plasmodium P. vivax

Circumsporozoite protein

Protein ID: -

Source:

Gene ID: CSP

No. of peptides: 157 peptides (Selected immunogenic peptides (15mers))

Amount: 25 μg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10⁸ cells)

Length: diverse

Sequence: see additional product documentation

Purity: Research Plus Grade: each peptide is guaranteed to be main product (HPLC/MS)

Appearance: white to off-white lyophilisate

Remark: The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.

Storage conditions: Store at - 20°C.

Application protocol: https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months

