

Product Data Sheet

Product:	PepMix Pan-Melanoma Select
Label name:	Pan-Melanoma Select
Product code:	PM-Pan-MELselect-2
Source:	Human Selected proteins
Protein ID:	n.a.
No. of peptides:	115 peptides (Defined HLA class I & II-restricted T-cell epitopes)
Amount:	15 nmol (approx. 25µg) / peptide for stimulation of 2,5 x 10 ⁸ cells
Length:	diverse
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months