

Product Data Sheet

Product:	PepMix Rabies virus (Glycoprotein)
Label name:	Rabies virus (Glycoprotein)
Product code:	PM-RABVH-GP
Source:	Rabies virus strain HEP-Flury Glycoprotein
Protein ID:	P19462
Gene ID:	G
No. of peptides:	129 peptides (15mers with 11aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	524 aa
Sequence:	MVPQVLLFAP LLVFPLCFGK FPIYTIPDKL GPWSPIDLHH LSCPNNLVVE DEGCTNLSGF SYMELKVGYI SAIKVNGFTC TGVVTEAETY TNFVGyvTTT FKRKHFRPTP DACRAAYNWK MAGDPRYEES LHNYPYDyHW LRTVKTTKES LVIISPSVTD LDPYDKSLHS RVFPGGNCSSG ITVSSTYCST NHDYTIWMPE NLRLGTSCDI FTHSRGKRAS KGDKTCGFVD ERGLYKSLKG ACKLKLCGVL GLRLMDGTWV AMQTSDETKW CPPGQLVNLH DFRSDEIEHL VEEELVKKRE ECLDALESIM TTKSVSFRRLL SHLRKLVPGF GKAYTIFNKT LMEADAHYKS VQTWNEIIPS KGCLRVGERC HPHVNGVFFN GIILGSDGHV LIPEMQSSLL QQHMELELESS VIPLMHPLAD PSTVFKDGDE VEDFVEVHLP DVHKQVSGVD LGLPKWGKYV LMIAGALIAL MLIIFLMTCC RRVNRPESTQ SNLGGTGRNV SVPSQSGKVI SSWESYKSGG ETRL
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C.
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/
Use recommended within 6 months	