

Product Data Sheet

Product:	PepMix Vaccinia virus (L3L)
Label name:	PepMix VACV (L3L)
Product code:	PM-VACV-L3L
Source:	Vaccinia Virus (VACV) strain Western Reserve Protein L3
Protein ID:	P07614
Gene ID:	VACWR090
No. of peptides:	85 peptides (15mers with 11 aa overlap)
Amount:	25 µg (approx. 15 nmol) / peptide (for stimulation of up to 2,5 x 10 ⁸ cells)
Length:	350 aa
Sequence:	MNTRTDVTND NIDKNPTKRG DKNIPGRNER FNDQNRFNND IPKPKPRLQP NQPPKQDNKC REENGDFINI RLCAYEKEYC NDGYLSPAYY MLKQVDDEEM SCWSELSSLV RSRKAVGFPL LKAAKRISHG SMLYFEQFKN SKVVRLTPQV KCLNDTVIFQ TVVILYSMYK RGIYSNEFCF DLVSIPRTNI VFSVNQLMFN ICTDILVVLV ICGNRLYRTN LPQSCYLNFI HGHETIARRG YEHSNYFFEW LIKNHISLLT KQTM DILKVK KYAIGAPVN RLLEPGTLVY VPKEDYYFIG ISLTDVSISD NRVVLFSTDG IVLEIEDFNI KHLFMAGEMF VRSQSSTIIV
Purity:	Development Grade: each peptide purified to > 70% (HPLC/MS)
Appearance:	white to off-white lyophilisate
Remark:	The peptides of this product are supplied as freeze-dried trifluoroacetate salts. For research use only.
Storage conditions:	Store at - 20°C
Application protocol:	https://jpt.com/support-contact/resources/application-protocols/

Use recommended within 6 months