

CD45R Mouse anti-Human Monoclonal (351C5) Antibody - LS-C134449 - LSBio

CatalogID:	LS-C134449
Target:	CD45R
Host	CD45R antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgM
Clone Name:	351C5
Immunogen Species:	CD45R antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45R antibody was raised against peripheral blood mononuclear cells.
Specificity:	CD45R isoforms of about 200-220 kd.
Reactivity:	Human
Purification:	Gel filtration
Presentation:	50 mM sodium phosphate, pH 7.5, 500 mM potassium chloride, 150 mM sodium chloride, 0.5 mg/ml Gentamicin sulfate
Recommended Storage:	Store at +4°C. DO NOT FREEZE.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 µg
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/25/2014

© 2014 LifeSpan BioSciences