

CD45RA Mouse anti-Human Monoclonal (PE) Antibody - LS-C21828 - LSBio

CatalogID:	LS-C21828
Target:	CD45RA
Host	CD45RA antibody was produced in Mouse
Clonality:	Monoclonal
Isotype:	IgG1
Conjugations:	Phycoerythrin (PE)
Immunogen Species:	CD45RA antibody was raised against Human
Antigen Type:	Cells
Immunogen:	CD45RA antibody was raised against human T lymphocytes.
Specificity:	Recognizes the high molecular weight form (220kD) of the human leucocyte common antigen (CD45RA). CD45RA is expressed by B lymphocytes, a subset of T lymphocytes, but not by granulocytes. Species cross-reactivity: Rhesus monkey. Does not cross-react with rat tissues.
Reactivity:	Human, Monkey
Purification:	Ion exchange chromatography
Presentation:	Lyophilized from PBS, pH 7.2, 0.1% sodium azide, 1% BSA.
Usage Summary:	Flow Cytometry: Neat.
Uses:	Flow Cytometry (Optimal dilution to be determined by the researcher)
Size:	100 tst
Requested From:	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 9/28/2014

© 2014 LifeSpan BioSciences