

**CD45R Mouse anti-Human Monoclonal Antibody - LS-C35620 - LSBio**

<b>CatalogID:</b>	LS-C35620
<b>Target:</b>	CD45R
<b>Host</b>	CD45R antibody was produced in Mouse
<b>Clonality:</b>	Monoclonal
<b>Isotype:</b>	IgM
<b>Immunogen Species:</b>	CD45R antibody was raised against Human
<b>Antigen Type:</b>	Cells
<b>Immunogen:</b>	CD45R antibody was raised against peripheral blood mononuclear cells.
<b>Specificity:</b>	Recognizes the CD45R isoforms of about 200-220kD.
<b>Reactivity:</b>	Human
<b>Purification:</b>	Protein A purified
<b>Presentation:</b>	PBS, pH 7.2, 5% glycerol, 0.2% BSA, 0.05% sodium azide.
<b>Usage Summary:</b>	Suitable for use in Flow Cytometry, Immunoprecipitation, Immunohistochemistry.
<b>Uses:</b>	IHC, Immunoprecipitation, Flow Cytometry (Optimal dilution to be determined by the researcher)
<b>Size:</b>	100 µg
<b>Requested From:</b>	Japan

Laboratory Reagent For In Vitro Research Use Only

Not for resale without prior written consent from LifeSpan BioSciences, Inc.

Created on 10/1/2014

© 2014 LifeSpan BioSciences