

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	01-25-703
Product-Name:	micromer ®
surface:	collagen
size:	7 µm
solid content:	25 mg/ml
quantity:	5 ml
Standard deviation:	< 5 % (C.V.)
Shape:	spherical
Density:	1.03 g/ccm
Stable in:	aqueous buffers
Not stable in:	organic solvents
Product form:	aqueous suspension
Particles per ml:	1.4*10E8
Particles per mg:	5.4*10E6
Additional remarks:	Storage at 2 - 8°C for 6 months