

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

| | |
|---------------------|---|
| Product-No.: | 30-02-503 |
| Product-Name: | micromer® -redF |
| surface: | COOH |
| size: | 5 µm |
| solid content: | 25 mg/ml |
| quantity: | 10 ml |
| Standard deviation: | < 5 % (C.V.) |
| Shape: | spherical |
| Density: | 1.03 g/ccm |
| Stable in: | aqueous buffers, methanol, ethanol, DMSO |
| Not stable in: | halogenated hydrocarbons, toluene, strong acidic solutions, eg 10% HCl |
| Product form: | aqueous suspension |
| Particles per ml: | 3.7*10E8 |
| Particles per mg: | 1.5*10E7 |
| Optical properties: | red fluorescent |
| Colour: | pink |
| Excitation: | 552 nm |
| Emission: | 580 nm |
| Additional remarks: | Storage at room temperature. |