

micromod Partikeltechnologie GmbH

Friedrich-Barnewitz-Str. 4
D-18119 Rostock
Germany

fon: +49 381 / 54 34 56 10
fax: +49 381 / 54 34 56 20
e-mail: info@micromod.de



Technical Datasheet

Product-No.:	40-01-701
Product-Name:	sicastar® -redF
surface:	NH₂
size:	70 nm
solid content:	25 mg/ml
quantity:	10 ml
Polydispersity index:	< 0.2
Shape:	spherical
Density:	2.0 g/ccm
Porosity:	nonporous
Stable in:	aqueous buffers, organic solvents
Not stable in:	hydrofluoric acid, strong basic media, e.g. 6 M NaOH
Product form:	aqueous suspension
Particles per ml:	7.0*10E13
Particles per mg:	2.8*10E12
Optical properties:	red fluorescent
Colour:	pink
Excitation:	569 nm
Emission:	585 nm
Additional remarks:	Storage at room temperature. Do not freeze!