

Product Datasheet

PPP2R5C Antibody (TQ11-1G6) [DyLight 680] NB100-501FR

Unit Size: 0.1 ml

Store at 4C in the dark.

www.novusbio.com



technical@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB100-501FR

Updated 10/1/2017 v.20.1

Earn rewards for product
reviews and publications.

Submit a publication at www.novusbio.com/publications

Submit a review at www.novusbio.com/reviews/destination/NB100-501FR



NB100-501FR

PPP2R5C Antibody (TQ11-1G6) [DyLight 680]

| Product Information | |
|---------------------|---|
| Unit Size | 0.1 ml |
| Concentration | Please see the vial label for concentration. If unlisted please contact technical services. |
| Storage | Store at 4C in the dark. |
| Clonality | Monoclonal |
| Clone | TQ11-1G6 |
| Preservative | 0.05% Sodium Azide |
| Isotype | IgG1 Kappa |
| Conjugate | DyLight 680 |
| Purity | Protein G purified |
| Buffer | 50mM Sodium Borate |

| Product Description | |
|-------------------------|---|
| Host | Mouse |
| Gene ID | 5527 |
| Gene Symbol | PPP2R5C |
| Species | Human, Mouse |
| Reactivity Notes | Human and mouse. |
| Specificity/Sensitivity | This is specific for PP2A (B56 gamma 3 subunit) protein. This does not recognize the B56 gamma 1 isoform of PP2A. |
| Immunogen | Human PPP2R5C [UniProt# Q13362] |
| Notes | Dylight (R) is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries. This conjugate is made on demand. Actual recovery may vary from the stated volume of this product. The volume will be greater than or equal to the unit size stated on the datasheet. |

| Product Application Details | |
|-----------------------------|--|
| Applications | Western Blot |
| Recommended Dilutions | Western Blot |
| Application Notes | Optimal dilution of this antibody should be experimentally determined. |





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NB100-501FR

| | |
|--------------------|---|
| NBP1-43319FR | Mouse IgG1 Kappa Light Chain Isotype Control (P3.6.2.8.1) [DyLight 680] |
| H00005527-P01-10ug | Recombinant Human PPP2R5C Protein |
| H00005520-P01-10ug | Recombinant Human PPP2R2A Protein |
| H00005527-Q01-10ug | Recombinant Human PPP2R5C Protein |

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our 100% guarantee, please visit www.novusbio.com/guarantee

Earn gift cards/discounts by submitting a review: www.novusbio.com/reviews/submit/NB100-501FR

Earn gift cards/discounts by submitting a publication using this product:
www.novusbio.com/publications

