

Product Datasheet

Histone H4 Antibody NB21-2088AF405

Unit Size: 50 Tests

Store at 4C in the dark.

www.novusbio.com



support@novusbio.com

Protocols, Publications, Related Products, Reviews, Research Tools and Images at:
www.novusbio.com/NB21-2088AF405

Updated 3/16/2014 v.20.1

NB21-2088AF405

Histone H4 [Monomethyl Lys20] Antibody [Alexa Fluor (R) 405]

Product Information	
Unit Size	50 Tests
Concentration	Please see the vial label for concentration. If unlisted please contact technical services.
Storage	Store at 4C in the dark.
Clonality	Polyclonal
Preservative	No Preservative
Conjugate	Alexa Fluor (R) 405
Purity	Immunogen affinity purified
Buffer	50 mM Sodium Borate

Product Description	
Host	Rabbit
Gene ID	121504
Gene Symbol	HIST4H4
Species	Human, Mouse, C. elegans
Species Reactivity	Human, mouse, and C. elegans. Predicted to react with many species including rat, chicken, Xenopus, Drosophila, and plant based on 100% sequence homology.
Marker	H4K20me1
Immunogen	Synthetic monomethylated peptide surrounding Lysine 20 of human Histone H4 [Swiss Prot P62805].
Notes	<p>Optimal dilution of this antibody should be experimentally determined.</p> <p>This product is provided under an intellectual property license from Life Technologies Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased product and components of the product only in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The sale of this product is expressly conditioned on the buyer not using the product or its components, or any materials made using the product or its components, in any activity to generate revenue, which may include, but is not limited to use of the product or its components: (i) in manufacturing; (ii) to provide a service, information, or data in return for payment; (iii) for therapeutic, diagnostic or prophylactic purposes; or (iv) for resale, regardless of whether they are resold for use in research. For information on purchasing a license to this product for purposes other than as described above, contact Life Technologies Corporation, 5791 Van Allen Way, Carlsbad, CA 92008 USA or outlicensing@lifetech.com.</p>

Product Application Details	
Applications	Western Blot, Chromatin Immunoprecipitation, Dot Blot, Immunocytochemistry/Immunofluorescence
Recommended Dilutions	Chromatin Immunoprecipitation, Dot Blot, Immunocytochemistry/Immunofluorescence, Western Blot



Application Notes

This Histone H4 K20me1 antibody is useful for Immunocytochemistry/Immunofluorescence ,ChIP, Dot blot, and Western blot where a band is seen ~13kDa.





Novus Biologicals USA

8100 Southpark Way, A-8
Littleton, CO 80120
USA
Phone: 303.730.1950
Toll Free: 1.888.506.6887
Fax: 303.730.1966
novus@novusbio.com

Novus Biologicals Canada

461 North Service Road West, Unit B37
Oakville, ON L6M 2V5
Canada
Phone: 905.827.6400
Toll Free: 855.668.8722
Fax: 905.827.6402
canada@novusbio.com

Novus Biologicals Europe

19 Barton Lane
Abingdon Science Park
Abingdon, OX14 3NB, United Kingdom
Phone: (44) (0) 1235 529449
Free Phone: 0800 37 34 15
Fax: (44) (0) 1235 533420
info@bio-techne.com

General Contact Information

www.novusbio.com
Technical Support: technical@novusbio.com
Orders: orders@novusbio.com
General: novus@novusbio.com

Products Related to NB21-2088AF405

NBP2-33376H	Blue Marker Antibody (6F4-F6) [HRP]
NB21-2088AF488	Histone H4 [Monomethyl Lys20] Antibody [Alexa Fluor (R) 488]
NB21-2103PEP	Histone H4 [Trimethyl Lys44] Peptide
NB21-2090PEP	Histone H4 [Trimethyl Lys20] Peptide

Limitations

This product is for research use only and is not approved for use in humans or in clinical diagnosis. Primary Antibodies are guaranteed for 1 year from date of receipt.

For more information on our guarantee, please visit www.novusbio.com/guarantee.

