

Mahogany Protein (1350-1376) (Mouse)

Catalog#	Standard size	Price	Order
003-86	100 µg	\$170.00	<input type="button" value="ORDER"/>

Product description

References

Catalog #	003-86
Standard Size	100 µg
Sequence	Thr - Val - Pro - Lys - Pro - Ile - Ala - Leu - Glu - Pro - Cys - Phe - Gly - Asn - Lys - Ala - Ala - Val - Leu - Ser - Val - Phe - Val - Arg - Leu - Pro - Arg
Species	Mouse
M.W	2923.6
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100µg of NET peptide.
HPLC Analysis	Purity= >80%