

Quality Control Record

Catalog No.: 007-22

Lot No.: 430432

Product: ELA, C-terminal / Prepro-ELA (44-54) (Human)

Sequence: Cys-Met-Pro-Leu-His-Ser-Arg-Val-Pro-Phe-Pro

Molecular Weight:
1283.57

Mass Spec Analysis:
Exhibits correct M.W.

HPLC Analysis:
Purity $\geq 95\%$

Appearance: White powder

Peptide Content:
This vial contains 100 μ g of NET peptide.

Solubility and Reconstitution:
Soluble in water.

Storage Conditions:
Store peptide in DRY form at 0-5°C.

Notes: For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions. Peptide contains free Cys and tends to dimerize.

Quality Control Manager: _____ **Date:** _____
Galina Citowsky