

NPQ-Gly (Human, Mouse, Bovine)

Catalog#	Standard size	Price	Order
020-10	200 µg	\$190.00	<input type="button" value="ORDER"/>

Product description

More information

[Spexin / NPQ \(Human, Mouse\)](#)

Catalog #	020-10
Standard Size	200 µg
Sequence	Asn - Trp - Thr - Pro - Gln - Ala - Met - Leu - Tyr - Leu - Lys - Gly - Ala - Gln - Gly
Species	Human, Mouse, Bovine
M.W	1677.94
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	The contents of this vial have been accurately determined to contain 200µg of NET peptide.
HPLC Analysis	Purity ≥95%

Thursday 20 December, 2012

11979858 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.