

## Spexin, Prepro (36-58) (Human)

Catalog#	Standard size	Price	Order
023-82	100 µg	\$190.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

[Spexin / NPO \(Human, Mouse\)](#)

<b>Catalog #</b>	023-82
<b>Family</b>	Spexin
<b>Standard Size</b>	100 µg
<b>Sequence</b>	Asn - Trp - Thr - Pro - Gln - Ala - Met - Leu - Tyr - Leu - Lys - Gly - Ala - Gln - Gly - Arg - Arg - Phe - Ile - Ser - Asp - Gln - Ser
<b>Species</b>	Human
<b>M.W</b>	2668.03
<b>Mass Spec Analysis</b>	Exhibits correct M.W.
<b>Solubility</b>	Soluble in water.
<b>Storage Condition</b>	Store peptide in DRY form at 0-5°C.
<b>Appearance</b>	White powder
<b>Content</b>	The contents of this vial have been accurately determined to contain 100 µg of NET peptide
<b>Additional articles</b>	For best results and reproducibility, rehydrate peptide just before use. Do not refreeze any unused portions. Do not attempt to weigh out smaller portions

Thursday 20 December, 2012

11980097 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.