

## Spexin, Prepro (73-116) (Human)

Catalog#	Standard size	Price	Order
023-83	100 µg	\$210.00	<input type="button" value="ORDER"/>

### Product description

### More information

### References

Spexin / NPO (Human, Mouse)

Catalog #	023-83
Family	Spexin
Standard Size	100 µg
Sequence	Ser - Pro - Asn - Pro - Gln - Leu - Leu - Thr - Ile - Pro - Glu - Ala - Ala - Thr - Ile - Leu - Leu - Ala - Ser - Leu - Gln - Lys - Ser - Pro - Glu - Asp - Glu - Glu - Lys - Asn - Phe - Asp - Gln - Thr - Arg - Phe - Leu - Glu - Asp - Ser - Leu - Leu - Asn - Trp
Species	Human
M.W	5013.59
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water.
Storage Condition	Store peptide in DRY form at 0-5°C.
Appearance	White powder
Content	The contents of this vial have been accurately determined to contain 100 µg of NET peptide
HPLC Analysis	Purity ≥95%
Additional articles	For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.

Thursday 20 December, 2012

11980099 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.