

Spexin / NPQ(53-70) (Human)

Catalog#	Standard size	Price	Order
023-84	100 µg	\$180.00	<input type="button" value="ORDER"/>

Product description

More information

[Spexin / NPQ \(Human, Mouse\)](#)

Catalog #	023-84
Family	Spexin
Standard Size	100 µg
Sequence	Phe - Ile - Ser - Asp - Gln - Ser - Arg - Arg - Lys - Asp - Leu - Ser - Asp - Arg - Pro - Leu - Pro - Glu
Species	Human
M.W	2159-38
Mass Spec Analysis	Exhibits correct M.W.
Solubility	Soluble in water
Storage Condition	Store peptide in DRY form at 0-5°C. For best results and reproducibility, rehydrate peptide just before use. Do not attempt to weigh out smaller portions. Do not refreeze any unused portions.
Appearance	White Powder
Content	100 µg NET peptide
HPLC Analysis	Purity >= 95%

Thursday 20 December, 2012

11980100 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.